

UPC-R80MD

A4 Self-Laminating Colour Printing Pack for UP-DR80MD



Overview

Self-laminating, A4 media for high durability and longevity.

The protective laminate coating is resistant to light, heat, water, humidity and fingerprints. Long-lasting prints for long-term storage.

One pack contains 2 ribbons and 2 rolls of media to yield 100 prints per pack. (4 packs per master carton)

Specifications

Specifications	
Type of Media	Self-Laminating Color Printing Pack
Contents	Paper roll for 50 prints x 2
Paper/Film Size	210 mm x 16 m 8 3/8 x 630 inches
Print Pages for One Roll	50 pages x 2/pack
Supported Printers	UP-DR80MD

Related products



UP-DR80MD

A4 Digital Colour Printer

Gallery

